

Product catalog

Dental Multi Unit Platforms Compatible



Zimmer 3.5 Platform Compatible

Straight multi unit abutment • D-type



1 mm
USMD-3701



2 mm
USMD-3702



3 mm
USMD-3703



4 mm
USMD-3704



5 mm
USMD-3705

Screw Driver Type: **Key for straight MUA D-type** •
Torque: **30 Ncm**

17° multi unit abutment • D-type



1 mm
UAMD-1701



2 mm
UAMD-1702



3 mm
UAMD-1703



4 mm
UAMD-1704

Supplied With: **Screw for angled MUA D-type Internal Hex RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



1 mm
UAMD-3001



2 mm
UAMD-3002



3 mm
UAMD-3003



4 mm
UAMD-3004

Supplied With: **Screw for angled MUA D-type Internal Hex RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type

UMSD-0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type



1 mm
UAMD-4501



2.5 mm
UAMD-4525

Supplied With: **Screw & Handle for angled MUA D-type Internal Hex RP**
Torque: **30 Ncm**

Zimmer 3.5 Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017


REF | UMTD-0001U



Transfer for MUA D-type


-  SUPPLIED WITH
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMSD-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMSD-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U


REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm


ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




Zimmer 3.5 Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment Internal Hex RP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment Internal Hex RP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment Internal Hex RP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment Internal Hex RP**

Zimmer 3.5 Platform Compatible

Screw for abutment

USWR-3783



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment

USWR-0002



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

USWR-T601



Screw Driver Type: **Torx 6**

Digital implant analog

UIAR-D001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

USBI-0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment Conical 11° RP**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	ULAR-0001C	ULAR-0002C	ULAR-0003C	ULAR-0004C	ULAR-0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-IHR0901C	LAR-IHR0902C	LAR-IHR0903C

15° DLOC attachment abutment


TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-IHR1501C	LAR-IHR1502C	LAR-IHR1503C

30° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-IHR3001C	LAR-IHR3002C	LAR-IHR3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Straumann BL RC 4.1-4.8 Compatible

Straight multi unit abutment • D-type



1 mm

SMD-SBR3701



2 mm

SMD-SBR3702



3 mm

SMD-SBR3703



4 mm

SMD-SBR3704



5 mm

SMD-SBR3705

Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



1 mm

AMD-SBR1701



2 mm

AMD-SBR1702



3 mm

AMD-SBR1703



4 mm

AMD-SBR1704

Supplied With: **Screw for angled MUA D-type Straumann®, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



1 mm

AMD-SBR3001



2 mm

AMD-SBR3002



3 mm

AMD-SBR3003



4 mm

AMD-SBR3004

Supplied With: **Screw for angled MUA D-type Straumann®, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type

MSD-SBR0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type



1 mm

AMD-SBR4501



2.5 mm

AMD-SBR4525

Supplied With: **Screw & Handle for angled MUA D-type Straumann®**
Torque: **30 Ncm**

Straumann BL RC 4.1-4.8 Compatible

Temporary sleeve for MUA D-type

-  **TORQUE**
15 Ncm
-  **SUPPLIED WITH**
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  **SUPPLIED WITH**
UMSD-0002U

REF | USCD-0001



Transfer for MUA D-type H15

-  **SUPPLIED WITH**
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  **SUPPLIED WITH**
UMSD-0017

REF | UMTD-0001U



Transfer for MUA D-type

-  **SUPPLIED WITH**
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMDS-0017



Screw TITLE

-  **TORQUE**
30 Ncm

REF | UMDS-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  **TORQUE**
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  **TORQUE**
30 Ncm
-  **SCREW DRIVER TYPE**
Hex, 1.27 mm

REF | UMDS-0002U



Screw for MUA sleeve D-type Torx6

-  **TORQUE**
30 Ncm
-  **SCREW DRIVER TYPE**
Hex, 1.27 mm

REF | UMDS-T603U



Digital analog MUA D-type

-  **SUPPLIED WITH**
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  **TORQUE**
15 Ncm
-  **SUPPLIED WITH**
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  **TORQUE**
15 Ncm
-  **SUPPLIED WITH**
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  **TORQUE**
30 Ncm
-  **SUPPLIED WITH**
UMSD-0002U


REF | USTD-0003C



Ratchet key for straight MUA D-type

-  **TORQUE**
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type

-  **TORQUE**
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730



Straumann BL RC 4.1-4.8 Compatible

Ti-Base UNO hex • H4



0.5 mm

TBU-SBR0400U



1 mm

TBU-SBR0401U



2 mm

TBU-SBR0402U



3 mm

TBU-SBR0403U



4 mm

TBU-SBR0404U

Supplied With: **Screw for abutment Straumann®**
Torque: **30 Ncm**

Ti-Base UNO hex • H6



0.5 mm

TBU-SBR0600U



1 mm

TBU-SBR0601U



2 mm

TBU-SBR0602U



3 mm

TBU-SBR0603U



4 mm

TBU-SBR0604U

Supplied With: **Screw for abutment Straumann®**
Torque: **30 Ncm**

Ti-Base UNO round • H4



0.5 mm

TBU-SBR1400U



1 mm

TBU-SBR1401U



2 mm

TBU-SBR1402U



3 mm

TBU-SBR1403U



4 mm

TBU-SBR1404U

Supplied With: **Screw for abutment Straumann®**
Torque: **30 Ncm**

Ti-Base UNO round • H6



0.5 mm

TBU-SBR1600U



1 mm

TBU-SBR1601U



2 mm

TBU-SBR1602U



3 mm

TBU-SBR1603U



4 mm

TBU-SBR1604U

Supplied With: **Screw for abutment Straumann®**
Torque: **30 Ncm**

Straumann BL RC 4.1-4.8 Compatible

Screw for abutment


SWR-SBR0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-SBRT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-SBR0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-SBR0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment Straumann®**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-SBR0001C	LAR-SBR0002C	LAR-SBR0003C	LAR-SBR0004C	LAR-SBR0005C

9° DLOC attachment abutment

TORQUE
30 Ncm

+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-SBR0901C	LAR-SBR0902C	LAR-SBR0903C

15° DLOC attachment abutment

TORQUE
30 Ncm

+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-SBR1501C	LAR-SBR1502C	LAR-SBR1503C

30° DLOC attachment abutment

TORQUE
30 Ncm

+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-SBR3001C	LAR-SBR3002C	LAR-SBR3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Nobel Active Regular Platform Compatible

Straight multi unit abutment • D-type




Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type NobelActive® RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type NobelActive® RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type

MSD-NAR0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type



Supplied With: **Screw & Handle for angled MUA D-type NobelActive® RP**
Torque: **30 Ncm**

Nobel Active Regular Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF

USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U

REF

USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF

UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF

UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF

UMTD-0002



Screw for MUA transfer D-type

REF

UMSD-0017



Screw TITLE

-  TORQUE
30 Ncm

REF

UMSD-2012



Analog MUA D-type

REF

UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF

UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF

UMSD-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF

UMSD-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF

UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF

USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF

USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF

USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM

H 6.3 mm

REF

UMDR-2716



Contra angle key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM

H 22 mm

REF

UMDC-2722


H 30 mm

UMDC-2730




Nobel Active Regular Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® RP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® RP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® RP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® RP**

Nobel Active Regular Platform Compatible

Screw for abutment

SWR-NAR0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-NART601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-NAR0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-NAR0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment NobelActive® RP**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-NAR0001C	LAR-NAR0002C	LAR-NAR0003C	LAR-NAR0004C	LAR-NAR0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NAR0901U	LAR-NAR0902U	LAR-NAR0903U

15° DLOC attachment abutment

TORQUE
30 Ncm
SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NAR1501U	LAR-NAR1502U	LAR-NAR1503U

30° DLOC attachment abutment

TORQUE
30 Ncm
SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NAR3001U	LAR-NAR3002U	LAR-NAR3003U

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Nobel Active Narrow Platform Compatible

Straight multi unit abutment • D-type



1 mm

2 mm

3 mm

4 mm

5 mm

Screw Driver Type: **Hex, 6.35 mm**

Torque: **30 Ncm**

SMD-NAN3701

SMD-NAN3702

SMD-NAN3703

SMD-NAN3704

SMD-NAN3705

17° multi unit abutment • D-type



1 mm

2 mm

3 mm

4 mm

AMD-NAN1701

AMD-NAN1702

AMD-NAN1703

AMD-NAN1704

Supplied With: **Screw for angled MUA D-type NobelActive® NP, Handle for angled MUA D-Type**

Torque: **30 Ncm**

30° multi unit abutment



1 mm

2 mm

3 mm

4 mm

AMD-NAN3001

AMD-NAN3002

AMD-NAN3003

AMD-NAN3004

Supplied With: **Screw for angled MUA D-type NobelActive® NP, Handle for angled MUA D-Type**

Torque: **30 Ncm**

Screw for angled MUA • D-type

MSD-NAN0007



Screw Driver Type: **Hex, 1.27 mm**

Nobel Active Narrow Platform Compatible

Temporary sleeve for MUA D-type

 TORQUE
15 Ncm

 SUPPLIED WITH
UMSD-0002U

REF

USTD-000



Castable sleeve for MUA D-type

 SUPPLIED WITH
UMSD-0002U

REF

USCD-0001




Transfer for MUA D-type H15

 SUPPLIED WITH
UMSD-0017

REF

UMTD-0002U



Transfer for MUA D-type H12

 SUPPLIED WITH
UMSD-0017

REF

UMTD-0001U



Transfer for MUA D-type

 SUPPLIED WITH
UMSD-2012

REF

UMTD-0002



Screw for MUA transfer D-type

REF

UMSD-0017



Screw TITLE

 TORQUE
30 Ncm

REF

UMSD-2012



Analog MUA D-type

REF

UMAD-0001U



Healing cap for MUA D-type

 TORQUE
15 Ncm

REF

UMHD-0004



Screw for MUA sleeve D-type U

 TORQUE
30 Ncm

 SCREW DRIVER TYPE
Hex, 1.27 mm

REF

UMSD-0002U



Screw for MUA sleeve D-type Torx6

 TORQUE
30 Ncm

 SCREW DRIVER TYPE
Hex, 1.27 mm

REF

UMSD-T603U



Digital analog MUA D-type

 SUPPLIED WITH
UIAR-S001

REF

UMAD-D001U



Scan body MUA level D-type

 TORQUE
15 Ncm

 SUPPLIED WITH
USBD-S001

REF

USBD-0001



Scan body MUA level D-type H12

 TORQUE
15 Ncm

 SUPPLIED WITH
USBD-S012

REF

USBD-0012



Sleeve for MUA D-type H4 C

 TORQUE
30 Ncm

 SUPPLIED WITH
UMSD-0002U

REF

USTD-0003C



Ratchet key for straight MUA D-type

 TORQUE
30 Ncm

ITEM

H 6.3 mm

REF

UMDR-2716



Contra angle key for straight MUA D-type

 TORQUE
30 Ncm

ITEM

H 22 mm

REF

UMDC-2722




H 30 mm

UMDC-2730




Nobel Active Narrow Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® NP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® NP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® NP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment NobelActive® NP**

Nobel Active Narrow Platform Compatible

Screw for abutment

SWR-NAN0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-NANT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-NAN0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-NAN0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment NobelActive® NP**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-NAN0001C	LAR-NAN0002C	LAR-NAN0003C	LAR-NAN0004C	LAR-NAN0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NAN0901U	LAR-NAN0902U	LAR-NAN0903U

15° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NAN1501U	LAR-NAN1502U	LAR-NAN1503U

30° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NAN3001U	LAR-NAN3002U	LAR-NAN3003U

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Megagen Anyridge Platform Compatible

Straight multi unit abutment • D-type



Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type Megagen Anyridge®, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type Megagen Anyridge®, Handle for angled MUA D-Type**
Torque: **30 Ncm**


Screw for angled MUA • D-type

MSD-MAR0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type



Supplied With: **Screw & Handle for angled MUA D-type Megagen Anyridge®**
Torque: **30 Ncm**

Megagen Anyrige Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF

USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U

REF

USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF

UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF

UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF

UMTD-0002



Screw for MUA transfer D-type

REF

UMSD-0017



Screw TITLE

-  TORQUE
30 Ncm

REF

UMSD-2012



Analog MUA D-type

REF

UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF

UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF

UMSD-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF

UMSD-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF

UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF

USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF

USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF

USTD-0003C



Ratchet key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM

H 6.3 mm

REF

UMDR-2716



Contra angle key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM

H 22 mm

REF

UMDC-2722




H 30 mm

UMDC-2730




Megagen Anyridge Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen Anyridge®**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen Anyridge®**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen Anyridge®**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen Anyridge®**

Megagen Anyridge Platform Compatible

Screw for abutment

SWR-MAR0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment

SWR-MART601

Torx6



Screw Driver Type: **Torx 6**

Extractor for angled

MED-MAR0001

D-type



Digital implant analog

IAD-MAR0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-MAR0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment Megagen Anyridge®**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-MAR0001C	LAR-MAR0002C	LAR-MAR0003C	LAR-MAR0004C	LAR-MAR0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-MAR0901U	LAR-MAR0902U	LAR-MAR0903U

15° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-MAR1501U	LAR-MAR1502U	LAR-MAR1503U

30° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-MAR3001U	LAR-MAR3002U	LAR-MAR3003U

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Megagen Anyone Platform Compatible

Straight multi unit abutment • D-type



Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type Megagen AnyOne®, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type Megagen AnyOne®, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type

MSD-MON0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type



Supplied With: **Screw & Handle for angled MUA D-type Megagen AnyOne®**
Torque: **30 Ncm**

Megagen Anyone Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMDS-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMDS-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U


REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




Megagen Anyone Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen AnyOne®**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen AnyOne®**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen AnyOne®**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment Megagen AnyOne®**

Megagen Anyone Platform Compatible

Screw for abutment

SWR-MON0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-MONT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-MON0001



Supplied With: **Nut for digital analog**
Screw Driver Type: **Hex, 1.27 mm**

Scan body implant level

SBI-MON0009




Supplied With: **Screw for abutment Megagen AnyOne®**
Torque: **5-10 Ncm**

**Download
LIBRARIES**



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-MON0001C	LAR-MON0002C	LAR-MON0003C	LAR-MON0004C	LAR-MON0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-MON0901C	LAR-MON0902C	LAR-MON0903C

15° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-MON1501C	LAR-MON1502C	LAR-MON1503C

30° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-MON3001C	LAR-MON3002C	LAR-MON3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003


Osstem | Hiossen Regular Platform Compatible

Straight multi unit abutment • D-type



Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type X11 RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type X11 RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type


MSD-UVR0007



Screw Driver Type: **Hex, 1.27 mm**



45° multi unit abutment • D-type



Supplied With: **Screw & Handle for angled MUA D-type X11 RP**
Torque: **30 Ncm**

Osstem | Hiossen Regular Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2002

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMSD-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMSD-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U


REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




Osstem | Hiossen Regular Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment X11 RP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment X11 RP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment X11 RP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment X11 RP**

Osstem | Hiossen Regular Platform Compatible

Screw for abutment

SWR-UVR0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-OSRT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-UVR0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-UVR0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment X11 RP**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-OSR0001C	LAR-OSR0002C	LAR-OSR0003C	LAR-OSR0004C	LAR-OSR0005C

9° DLOC attachment abutment

TORQUE
30 Ncm

SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-OSR0901C	LAR-OSR0902C	LAR-OSR0903C

15° DLOC attachment abutment

TORQUE
30 Ncm

SUPPLIED WITH
USWR-3783




ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-OSR1501C	LAR-OSR1501C	LAR-OSR1503C

30° DLOC attachment abutment

TORQUE
30 Ncm

SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-OSR3001C	LAR-OSR3002C	LAR-OSR3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Osstem | Hiossen Mini Platform Compatible

Straight multi unit abutment • D-type



Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type Conical 11° MP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type Conical 11° MP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type

MSD-OSM0007



Screw Driver Type: **Hex, 1.27 mm**

Osstem | Hiossen Mini Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMDS-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMDS-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




Osstem | Hiossen Mini Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment Conical 11° MP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment Conical 11° MP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment Conical 11° MP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment Conical 11° MP**

Osstem | Hiossen Mini Platform Compatible

Screw for abutment

UOSM-0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-OSMT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-OSM0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-OSM0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment Conical 11° MP**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-OSM0001C	LAR-OSM0002C	LAR-OSM0003C	LAR-OSM0004C	LAR-OSM0005C

9° DLOC attachment abutment


TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-OSM0901C	LAR-OSM0902C	LAR-OSM0903C

15° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-OSM1501C	LAR-OSM1502C	LAR-OSM1503C

30° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-OSM3001C	LAR-OSM3002C	LAR-OSM3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Neodent GM Platform Compatible


Straight multi unit abutment • D-type



1 mm SMD-NER3701 2 mm SMD-NER3702 3 mm SMD-NER3703 4 mm SMD-NER3704 5 mm 5129.1005

Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



1 mm AMD-NER1701 2 mm AMD-NER1702 3 mm AMD-NER1703 4 mm AMD-NER1704

Supplied With: **Screw for angled MUA D-type Neodent GM®, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



1 mm AMD-NER3001 2 mm AMD-NER3002 3 mm AMD-NER3003 4 mm AMD-NER3004

Supplied With: **Screw for angled MUA D-type Neodent GM®, Handle for angled MUA D-Type**
Torque: **30 Ncm**


Screw for angled MUA • D-type

MSD-NER0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type



1 mm AMD-NER4501 2.5 mm AMD-NER4525

Supplied With: **Screw & Handle for angled MUA D-type Neodent GM®**
Torque: **30 Ncm**

Neodent GM Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMDS-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMDS-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U


REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




Neodent GM Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment Neodent GM®**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment Neodent GM®**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment Neodent GM®**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment Neodent GM®**

Neodent GM Platform Compatible

Screw for abutment

SWR-NER0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-NERT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-NER0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-NER0009




Torque: **5-10 Ncm**
Supplied With: **Screw for abutment Neodent GM®**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-NER0001C	LAR-NER0002C	LAR-NER0003C	LAR-NER0004C	LAR-NER0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NER0901C	LAR-NER0902C	LAR-NER0903C

15° DLOC attachment abutment


TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NER1501C	LAR-NER1502C	LAR-NER1503C

30° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-NER3001C	LAR-NER3002C	LAR-NER3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

Dentium Platform Compatible

Straight multi unit abutment • D-type



Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type Dentium® RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type Dentium® RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type

MSD-DTR0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type



Supplied With: **Screw & Handle for angled MUA D-type Dentium® RP**
Torque: **30 Ncm**

Dentium Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMSD-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMSD-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U


REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




Dentium Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment Dentium® RP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment Dentium® RP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment Dentium® RP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment Dentium® RP**

Dentium Platform Compatible

Screw for abutment

SWR-DTR0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-DTRT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-DTR0001



Screw Driver Type: **Hex, 1.27 mm** •
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-DTR0009




Torque: **5-10 Ncm** •
Supplied With: **Screw for abutment Dentium® RP**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-DTR0001C	LAR-DTR0002C	LAR-DTR0003C	LAR-DTR0004C	LAR-DTR0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-DTR0901C	LAR-DTR0902C	LAR-DTR0903C

15° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-DTR1501C	LAR-DTR1502C	LAR-DTR1503C

30° DLOC attachment abutment

TORQUE
30 Ncm
 SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-DTR3001C	LAR-DTR3002C	LAR-DTR3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

MIS C1/V3 SP Platform Compatible

Straight multi unit abutment • D-type



Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type MIS® V3 RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type MIS® V3 RP, Handle for angled MUA D-Type**
Torque: **30 Ncm**


Screw for angled MUA • D-type

MSD-M3R0007



Screw Driver Type: **Hex, 1.27 mm**

45° multi unit abutment • D-type




Supplied With: **Screw & Handle for angled MUA D-type MIS® V3 RP**
Torque: **30 Ncm**

MIS C1/V3 SP Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMDS-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMDS-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMDS-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U


REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




MIS C1/V3 SP Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 RP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 RP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 RP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 RP**

MIS C1/V3 SP Platform Compatible

Screw for abutment

SWR-M3R0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-M3RT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-M3R0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-M3R0009




Torque: **5-10 Ncm** •
Supplied With: **Screw for abutment MIS® V3 RP**

Download LIBRARIES



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-M3R0001C	LAR-M3R0002C	LAR-M3R0003C	LAR-M3R0004C	LAR-M3R0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-M3R0901C	LAR-M3R0902C	LAR-M3R0903C

15° DLOC attachment abutment


TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-M3R1501C	LAR-M3R1502C	LAR-M3R1503C

30° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-M3R3001C	LAR-M3R3002C	LAR-M3R3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003

MIS V3 NP Platform Compatible

Straight multi unit abutment • D-type



Screw Driver Type: **Key for straight MUA D-type**
Torque: **30 Ncm**

17° multi unit abutment • D-type



Supplied With: **Screw for angled MUA D-type MIS® V3 NP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

30° multi unit abutment



Supplied With: **Screw for angled MUA D-type MIS® V3 NP, Handle for angled MUA D-Type**
Torque: **30 Ncm**

Screw for angled MUA • D-type

MSD-M3N0007



Screw Driver Type: **Hex, 1.27 mm**

MIS V3 NP Platform Compatible

Temporary sleeve for MUA D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-000



Castable sleeve for MUA D-type

-  SUPPLIED WITH
UMSD-0002U


REF | USCD-0001



Transfer for MUA D-type H15

-  SUPPLIED WITH
UMSD-0017

REF | UMTD-0002U



Transfer for MUA D-type H12

-  SUPPLIED WITH
UMSD-0017


REF | UMTD-0001U



Transfer for MUA D-type

-  SUPPLIED WITH
UMSD-2012

REF | UMTD-0002



Screw for MUA transfer D-type

REF | UMSD-0017



Screw TITLE

-  TORQUE
30 Ncm

REF | UMSD-2012



Analog MUA D-type

REF | UMAD-0001U



Healing cap for MUA D-type

-  TORQUE
15 Ncm

REF | UMHD-0004



Screw for MUA sleeve D-type U

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-0002U



Screw for MUA sleeve D-type Torx6

-  TORQUE
30 Ncm
-  SCREW DRIVER TYPE
Hex, 1.27 mm

REF | UMSD-T603U



Digital analog MUA D-type

-  SUPPLIED WITH
UIAR-S001

REF | UMAD-D001U



Scan body MUA level D-type

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S001

REF | USBD-0001



Scan body MUA level D-type H12

-  TORQUE
15 Ncm
-  SUPPLIED WITH
USBD-S012

REF | USBD-0012



Sleeve for MUA D-type H4 C

-  TORQUE
30 Ncm
-  SUPPLIED WITH
UMSD-0002U

REF | USTD-0003C



Ratchet key for straight MUA D-type

-  TORQUE
30 Ncm

ITEM	H 6.3 mm
REF	UMDR-2716



Contra angle key for straight MUA D-type


-  TORQUE
30 Ncm

ITEM	H 22 mm	H 30 mm
REF	UMDC-2722	UMDC-2730




MIS V3 NP Platform Compatible

Ti-Base UNO hex • H4




Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 NP**

Ti-Base UNO hex • H6




Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 NP**

Ti-Base UNO round • H4



Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 NP**

Ti-Base UNO round • H6



Torque: **30 Ncm** • Supplied With: **Screw for abutment MIS® V3 NP**

MIS V3 NP Platform Compatible

Screw for abutment

SWR-M3N0001



Screw Driver Type: **Hex, 1.27 mm**

Screw for abutment • Torx6

SWR-M3NT601



Screw Driver Type: **Torx 6**

Digital implant analog

IAD-M3N0001



Screw Driver Type: **Hex, 1.27 mm**
Supplied With: **Nut for digital analog**

Scan body implant level

SBI-M3N0009




Torque: **5-10 Ncm** •
Supplied With: **Screw for abutment MIS® V3 NP**

**Download
LIBRARIES**



DLOC attachment abutment

SCREW DRIVER TYPE
Hex, 1.27 mm



ITEM	GH 1 mm	GH 2 mm	GH 3 mm	GH 4 mm	GH 5 mm
REF	LAR-M3N0001C	LAR-M3N0002C	LAR-M3N0003C	LAR-M3N0004C	LAR-M3N0005C

9° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-M3N0901C	LAR-M3N0902C	LAR-M3N0903C

15° DLOC attachment abutment


TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-M3N1501C	LAR-M3N1502C	LAR-M3N1503C

30° DLOC attachment abutment

TORQUE
30 Ncm
+ SUPPLIED WITH
USWR-3783



ITEM	GH 1 mm	GH 2 mm	GH 3 mm
REF	LAR-M3N3001C	LAR-M3N3002C	LAR-M3N3003C

Housing cap for DLOC



REF | ULMR-0001

Silicone cap for DLOC lab



REF | ULCR-0004

Silicone cap for DLOC standard



REF | ULCR-0002

Silicone cap for DLOC soft



REF | ULCR-0001

Silicone cap for DLOC strong



REF | ULCR-0003



Smart Implants

13 Aluf Yacov Perry Str.
Office 304,
Rehovot 7639302
Israel

Info@implant-smart.com
www.Implant-smart.com

 +972(0)54-3575050